

TITLEVIEWER

CAM\_TITLE

:  
OnEvent( "TITLEVIEWER" , "\_id\_" , "\_\_" )

:  
DoReact( "CAM\_TITLE" , "\_id\_" , "\_\_" );

:  
4 , .

OnEvent( "TITLEVIEWER" , "1" , "GO\_VIDEO" )  
{  
    DoReact( "MONITOR" , "4" , "ARCH\_FRAME\_TIME" , "cam<" +cam+ ">,date<" +date+ ">,time<" +time+ ">" );  
    DoReact( "MONITOR" , "4" , "KEY\_PRESSED" , "key<PLAY>" );  
}

:  
1.

OnEvent( "MACRO" , "1" , "RUN" )  
{  
    DoReact( "CAM\_TITLE" , "2" , "REINDEX" );  
}